# Dr. STORZ® Aloe Vera 97,5 % Gel\*

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### 1 Anwendungszweck

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## Instructions for Use

### 1 Therapeutic purpose

*Dr. Storz® Aloe Vera 97.5 % Gel\** is a medical device and applied to provide relief from itching.

The product makes use of the physical evaporation coolness of the ingredients to cool itching skin regions. The cooling reduces itching and thus stops the desire to scratch. This supports the symptomatic therapy of the causative disorder.

### 2 Therapeutic indications

*Dr. Storz® Aloe Vera 97.5 % Gel\** is ready-to-use, is intended for topical application, or regional application to itching skin, caused by insect stings or 1<sup>st</sup> degree sunburn. The product is suitable for application as required or for constant application at home and for clinical use for the skin of children from 3 years, adolescents and adults.

### 3 Undesirable effects

In rare cases the product can elicit allergic skin reactions.

### 4 Histocompatibility and interactions

There is no data indicating a histoincompatibility to *Dr. Storz® Aloe Vera 97.5 % Gel\**. *Dr. Storz® Aloe Vera 97.5 % Gel\** should not be applied simultaneously and on the same area with other creams, gels and ointments.

### 5 Contraindications

*Dr. Storz® Aloe Vera 97.5 % Gel\** may not be applied to skin diseases, open wounds, sunburn more severe than 1<sup>st</sup> degree nor when hypersensitivity to one of the ingredients is known.

#### Pregnancy and lactation

During breastfeeding, the use of *Dr. Storz® Aloe Vera 97.5 % Gel\** on or near the nipples should be avoided.

### 6 General special warnings and special precautions for use

*Dr. Storz® Aloe Vera 97.5 % Gel\** is exclusively intended for external application and may not enter the body by ingestion, injection or infusion.

*Dr. Storz® Aloe Vera 97.5 % Gel\** should not get into the eyes. In case of eye contact, please rinse thoroughly with water and consult an ophthalmologist if symptoms persist.

#### Attention:

Keep this medical device out of the reach and sight of children. Store in a safe place!

#### Attention:

Before taking *Dr. Storz® Aloe Vera 97.5 % Gel\** from the package, check the latter for integrity.

### 7 Product properties

Insect stings and sunburn often cause tormenting itching. The same applies to dry, aging skin which is often accompanied by intensive itching.

The physical effect of cooling can reduce itching of the skin. Particularly moisture saturated preparations provide relief from itching with their evaporation coolness. The ingredient Aloe Vera has a high hydropexis quality and additionally enhances cooling.

*Dr. Storz® Aloe Vera 97.5 % Gel\** was tested on 50 subjects in a double-blind clinical study at the Charité Berlin in connection for its cooling properties in connection with insect stings. A significant increase of the cooling effect was found.<sup>1</sup>

<sup>1</sup> Randomized, double-blind, partly open, controlled clinical trial on the efficacy and compatibility of dioxopromethazin hydrochloride (Prothanon® Gel), camphor (Vaopin® cooling gel) and Aloe Vera (*Dr. Storz® Aloe Vera 97.5 % Gel\**) in topical application compared to a placebo and / or non-treatment of test persons with provoked insect stings (Aedes aegyptii), 2004, Charité Berlin, Division of Evidence Based Medicine, Prof. Dr. med. B. Rzany)

### 8 Product application

*Dr. Storz® Aloe Vera 97.5 % Gel\** is applied as a covering layer and rubbed in thoroughly to the skin regions to be treated if and when required.

If the product is kept in the fridge (not below 5 °C), its cooling effect is enhanced. Close opened containers tightly after use. Used up containers and packaging are disposed of in domestic waste, preferably separated.

### 9 Physical product data

#### Composition:

Aloe barbadensis (concentrate reconstituted in cleaned water), Glycerol 85 %, Phenoxyethanol, Carbopol® ETD 2020, sodium edetate, sodium hydroxide, Tagat L2, perfume oil Chevrefeuille

#### Appearance, Odour:

Colourless to slightly yellow gel of refreshing bloomy odour.

#### Shelf-life:

30 months

The minimum shelf-life is stamped onto the bottom seam of the tube.

#### Storage conditions:

The product is stored in the intact packaging, dry at a storage temperature of 5 °C to 25 °C.

### 10 Supply specifications

3 g Sachet

100 g tube, packed in folding box

200 g tube, packed in folding box

500 g dispenser bottle

### 11 Manufacturer / Distributor

Manufacturer and distributor of the medical device:

 Esteve Pharmaceuticals GmbH  
Hohenzollerndamm 150-151  
14199 Berlin  
Germany  
Tel.: 030 338427-0  
Fax: 030 897480-45  
E-mail info.germany@esteve.com

Manufacturer

Temperature limits

Shelf life

Batch definition/Batch number

Note Instructions for Use

#### Note:

Please read instructions for use!

This leaflet was last approved in 06/2022

\* This product is manufactured from an Aloe Vera concentrate.

CE 0483

114118628

ESTEVE